

Comparison Chart of Equine Athletic Recovery Modalities

	EQ Press	Game Ready	Hidez Compression Suit	Bemer Horse Set	Cold Salt Water Spa
Dynamic (directional) pneumatic compression	✓	—*	—**	—	—
Patented smart technology to prevent backflow***	✓	—	—	—	—
Multiple chambers inflating providing a milking effect from bottom up the legs	✓	—	—	—	—
Precisely controlled pressure and treatment times	✓	—	—	—	—
Full Limb (Foot to elbow or stifle)	✓	—	✓	—	✓
Able to treat 1,2 or all 4 limbs at same time	✓	✓	—	—	—
All batteries/pumps incorporated (stand alone device)	✓	—	N/A	✓	N/A

*Non-dynamic pneumatic compression

**Static/passive compression

***A pressure gradient is created from the bottom chamber towards the top to promote forward flow against gravity and eliminate any possibility of backflow.

Note: Modalities such as shock wave, laser, and therapeutic ultrasound are intended for the focal treatment of a specific area of the body or a specific injury and therefore are not listed in the table above.